




UNZILA INDUSTRY

**ALL TYPES
OF DENTAL
INSTRUMENTS**



contact us

Mobile No : 0092 3086403561


Web Site : www.unzilaindustry.com

E-mail : unzilaindustry@gmail.com

Address : Ugoki road near sobedar street Sialkot




TABLE OF CONTENTS



01

Extracting Forceps


Page 01-31



02

Root Elevators


Page 32-50



03

Filling Instruments

Page 51-55



04

Dental Probes


Page 56-59



05

Periodontia Instruments
Mallets

Page 60-67



06


Curettes
Scalers

Page 68-74

TABLE OF CONTENTS

Page 75-77 Modelling Instruments
Measuring Instruments

07




Page 78-79 Wax / Plaster Knives

08




Page 80-85 Diagnostic (Explorers)
Endodontic Instruments

09



Page 86-87 Excavators

10



Page 88-89 Handles For Mouth Mirrors
Mouth Mirrors

11




Page 90-93 Spatulas

12




TABLE OF CONTENTS



13

Needle Holders &
Saliva Ejectors

Page 94-100



14

Syringes & Scissors


Page 101-109



15

Dental Tweezers


Page 110-114



16

Rubber Dam Clamps

Page 115-117



17

Matrix Retainers &
Matrix Bands

Page 118-119









18

Crown Removers, Scissors,
Forceps, Amalgam Guns & Carriers

Page 120-122

TABLE OF CONTENTS

Page 123	Copper Ring Removing Pliers & Paper Articulating Forceps	19	
Page 124-127	Rongeurs	20	
Page 128-129	Bone Elevators, Files & Chisels	21	
Page 130-151	Pliers For Orthodontics, Mouth Gags, Tongue Depressors Lip and Cheek Retractors and Implantology Instrumetns	22	
Page 152-155	Impression Trays	23	
Page 156-157	Articulators	24	
Page 158-160	Hollowares	25	

Page 1-31

- Extracting Forceps - English Pattern



UI 1-5101
Fig 8
lower premolars



UI 1-5102
Fig 13
lower premolars



UI 1-5103
Fig 17
upper molars, right



UI 1-5104
Fig 18
upper molars, left




UI 1-5105
Fig 18 A
upper molars






UI-1-5106
Fig 19
upper wisdoms



UI-1-5107
Fig 20
lower wisdoms



UI-1-5108
Fig 21
lower molars



UI-1-5109
Fig 22
lower molars



UI-1-5110
Fig 22 A
Albrecht
lower molars



UI-1-5111
Fig 33 A
lower roots



UI-1-5112
Fig 33 B
lower roots



UI-1-5113
Fig 33 L
lower roots




UI-1-5114
Fig 33 L
lower roots




UI-1-5115
Fig 33 S
lower roots, for children




UI-1-5116
Fig 44
upper roots




UI-1-5117
Fig 45
lower roots



UI-1-5118
Fig 12
upper roots




UI-1-5119
Fig 51 A
upper roots




UI-1-5120
Fig 51 L
upper roots




Relax
UI-1-5121
Fig 51
upper roots




Relax
UI-1-5122
Fig 51 A
upper roots



Relax
UI-1-5123
Fig 51 L
upper roots



Relax
UI-1-5124
Fig 67
upper wisdoms



Relax
UI-1-5125
Fig 67 A
upper wisdoms



UI-1-5126
Fig 37
Upper incisors



UI-1-5127
Fig 39
Upper molars



UI-1-5128
Fig 39 R
Upper molars, right



UI-1-5129
Fig 39 L
Upper molars, left



UI-1-5130
Fig 515
Upper roots




Klein
UI-1-5132
Fig 137
upper incisors



Klein
UI-1-5131
Fig 6
lower molars



Klein
UI-1-5133
Fig 5
lower incisors



Klein
UI-1-5134
Fig 139
uper premolars




UI-1-5135
Fig 43 A
Fractured or decayed
upper teeth



UI-1-5136
Fig 44
Upper roots



UI-1-5137
Fig 45
Lower roots



UI-1-5138
Fig 46
Lower roots



UI-1-5139
Fig 46 L
Lower roots



UI-1-5140
Fig 48
Upper roots



WITH SERRATED
TIPS



UI-1-5141
Fig 51
Upper roots



WITH SERRATED
TIPS



UI-1-5142
Fig 49
Upper roots



WITH SERRATED
TIPS



UI-1-5143
Fig 51 L
Upper roots



WITH SERRATED
TIPS



UI-1-5144
Fig 51 A
Upper roots



WITH SERRATED
TIPS



UI-1-5145
Fig 51 LX
Upper roots



1046
Fig 52
Upper roots and
canines



UI-1-5147
Fig 54
For separating
upper molars



UI-1-5148
Fig 55
For separating
upper molars




UI-1-5149
Fig 56
For separating
lower molars





UI-1-5155
Fig 66 L
Upper molars, left



UI-1-5156
Fig 66 R
Upper molars, right



UI-1-5157
Fig 67
Upper third molars




UI-1-5158
Fig 67 A
Upper third molars



UI-1-5159
Fig 68
Lower roots

Dental Instruments Vector Set



Excavator

Tweezers

Syringe

Dental
Nippers

Mouth
Mirror

Excavator



UI-1-5160
Fig 67½ L
Upper third molars, left



UI-1-5161
Fig 67½ R
Upper third molars, right



UI-1-5162
Fig 73
Lower molars



UI-1-5163
Fig 74
Lower roots



UI-1-5164
Fig 74 M
Lower roots



UI-1-515
Fig 74 N
Lower roots




WITH SERRATED
TIPS

UI-1-5166
Fig 75
Lower canines and
premolars




WITH SERRATED
TIPS

UI-1-5167
Fig 79
Lower third molars



WITH SERRATED
TIPS

UI-1-5168
Fig 79 A
Lower third molars



WITH SERRATED
TIPS

UI-1-5169
Fig 79 C
Lower third molars



UI-1-5170
Fig 86 W
Lower third molars



UI-1-5171
Fig 86 1/2 L
Lower third molars, left



UI-1-5172
Fig 86 1/2 R
Lower third molars, right



UI-1-5173
Fig 86 C
Lower molars



UI-1-5174
Fig 87
Lower decayed molars
or broken-down crowns




WITH SERRATED
TIPS

UI-1-5175
Fig 89
Upper molars, right



WITH SERRATED
TIPS

UI-1-5176
Fig 90
Upper molars, left



WITH SERRATED
TIPS

UI-1-5177
Fig 99
Broken-down lower molars



WITH SERRATED
TIPS

UI-1-5178
Fig 107
Upper canines



UI-1-5179
Fig 121
Upper third molars



UI-1-5180
Fig 145
Lower roots



WITH SERRATED
TIPS



UI-1-5181
Fig 151½
Upper roots



WITH SERRATED
TIPS



UI-1-5182
Fig 166
Universal extracting forceps for
upper teeth, in particular molars
and premolars



WITH SERRATED
TIPS



UI-1-5183
Fig 167
Universal extracting forceps for
lower teeth, in particular molars
and premolars




WITH SERRATED
TIPS




UI-1-5184
Fig 168
Upper canines



UI-1-5185
Fig 169
Lower canines



UI-1-5186
Fig 346
Lower roots



UI-1-5187
Fig 349
Upper roots




UI-1-5188
Fig 359
Lower roots



UI-1-5189
Fig 1
Universal extracting forcep for
upper incisors and canines



UI-1-5190
Fig 2
Upper molars



UI-1-5191
Fig 3
Lower incisors, premolar
and roots

WITH SERRATED
TIPS



UI-1-5192
Fig 4
Lower molars






UI-1-5193
Fig 1
Upper incisors
and canines



UI-1-5194
Fig 2
Upper incisors
and canines



UI-1-5195
Fig 4
Lower incisors
and canines



UI-1-5196
Fig 7
Upper premolars



UI-1-5197
Fig 13
Lower premolars




UI-1-5198
Fig 17
Upper molars, right



WITH SERRATED
TIPS

UI-1-5199
Fig 18
Upper molars, left




WITH SERRATED
TIPS

UI-1-5200
Fig 21
Lower molars



WITH SERRATED
TIPS

UI-1-5201
Fig 22
Lower molars



WITH SERRATED
TIPS

UI-1-5202
Fig 33
Lower roots



UI-1-5203
Fig 45
Lower roots



UI-1-5204
Fig 51 A
Upper roots



UI-1-5205
Fig 67 A
Upper third molars



UI-1-5206
Fig 74
Lower roots



UI-1-5207
Fig 79
Lower third molars



UI-1-5208
Fig 86 C
Lower molars



UI-1-5209
Fig 359
Lower roots



UI-1-5210
Fig 89
Upper molars, right



WITH SERRATED
TIPS



UI-1-5211
Fig 90
Upper molars, left



WITH SERRATED
TIPS



UI-1-5212
Fig 145
Lower roots



WITH SERRATED
TIPS



UI-1-5213
Fig 151 1/2
Upper roots



WITH SERRATED
TIPS



UI-1-5214
Fig 349
Upper roots



WITH SERRATED
TIPS



UI-1-5215
Fig 5
Lower third molars



UI-1-5216
Fig 6
Separating forcep



UI-1-5217
Fig 105
Upper molars



UI-1-5218
Fig 16
Lower molars



UI-1-5219
Fig 17
Lower molars



UI-1-5220
Fig 18 L
Upper molars, left



UI-1-5221
Fig 18 R
Upper molars, right



UI-1-5222
Fig 23
Lower molars



UI-1-5223
Fig 53 R
Upper molars, right



UI-1-5224
Fig 53 L
Upper molars, left



UI-1-5225
Fig 65
Upper incisors and roots



UI-1-5226
Fig 69
Upper and lower roots



UI-1-5227
Fig 88 R
Upper molars, right



UI-1-5228
Fig 88 L
Upper molars, left



UI-1-5229
Fig 101
Deciduous teeth



UI-1-5230
Fig 150
Upper incisors, premolars, roots



UI-1-5231
Fig 150 S
Upper teeth (for children)



UI-1-5232
Fig 151
Upper incisors, premolars, roots



UI-1-5233
Fig 151S
Lower teeth (for children)



UI-1-5234
Fig 222
Lower third molars




UI-1-5235
Fig 1
upper incisors and bicuspid



UI-1-5236
Fig 1 B
upper roots, anterior



Physick
UI-1-5237
Fig 5
lower wisdoms



UI-1-5238
Fig 6
lower molars, splitting forceps



UI-1-5239
Fig 10 H
upper molars



UI-1-5240
Fig 135
Lower premolars



WITH SERRATED
TIPS



UI-1-5241
Fig 225
Lower molars



UI-1-5242
Fig 305
Upper incisors
and roots



UI-1-5243
Fig 33 S
Lower roots



WITH SERRATED
TIPS



UI-1-5244
Fig 38
Lower incisors
and canines



WITH SERRATED
TIPS



UI-1-5245
Fig 1
Upper incisors



UI-1-5246
Fig 2
Upper premolars



UI-1-5247
Fig 3
Upper molars



UI-1-5248
Fig 4
Upper roots



UI-1-5249
Fig 5
Lower incisors



UI-1-5250
Fig 6
Lower molars



UI-1-5251
Fig 7
Lower roots



UI-1-5252
Fig 1
Lower roots



UI-1-5253
Fig 2
Upper roots



UI-1-5255
Fig 2
Roots fragments



UI-1-5254
Fig 1
Roots fragments


Page 32-50

- Root Elevators



lower jaw, left side

lower jaw, right side



Root elevator, spiral-shaped, for extremely facilitated extraction of the molars in the lower jaw. The slender, conical points allow deep penetration between tooth and gum. Owing to the hook-shaped working parts, slipping of the elevator is avoided.



UI-2-5101



UI-2-5102



Bein



Bein
UI-2-5103
Fig 1



Bein
UI-2-5104
Fig 2



Bein
UI-2-5105
Fig 0



Bein
UI-2-5106
Fig 3



Bein
UI-2-5107
Fig 4



Bein
UI-2-5108
Fig 5



Bein
UI-2-5109
Fig 6



Seldin



Seldin
UI-2-5110
Fig 1L



Seldin
UI-2-5111
Fig 1R



Seldin
UI-2-5112
Fig 2



Seldin
UI-2-5113
Fig 3



Seldin
UI-2-5114
Fig 4L



Seldin
UI-2-5115
Fig 4R



**Cryer
UI-2-5116
Fig 3**



**Cryer
UI-2-5117
Fig 4**



**Cryer
UI-2-5118
Fig 12**



**Cryer
UI-2-5119
Fig 27**



**Cryer
UI-2-5120
Fig 28**



**Cryer
UI-2-5121
Fig 27 S**



**Cryer
UI-2-5122
Fig 28 S**



Apical



Apical
UI-2-5123
Fig 1





Apical
UI-2-5124
Fig 2



Apical
UI-2-5125
Fig 3



Apical
UI-2-5126
Fig 4





UI-2-5132
Fig 1



UI-2-5133
Fig 2



UI-2-5134
Fig 3



Flohr
UI-2-5135
Fig 1



Flohr
UI-2-5136
Fig 2



Flohr
UI-2-5137
Fig 3



Kopp
UI-2-5138
Fig 1



Kopp
UI-2-5139
Fig 2



Kopp
UI-2-5144
Fig 1



Kopp
UI-2-5145
Fig 2



Kopp



Kopp
UI-2-5140
Fig 3



Kopp
UI-2-5141
Fig 4



Kopp



Kopp
UI-2-5142
Fig 5



Kopp
UI-2-5143
Fig 6




Kopp
UI-2-5146
Fig 3




Kopp
UI-2-5147
Fig 4



Hylin
UI-2-5148
Fig 1



Bein
UI-2-5150
Fig 1



Bein
UI-2-5151
Fig 2




Hylin



Hylin
UI-2-5149
Fig 2



Bein



Bein
UI-2-5152
Fig 3



Bein
UI-2-5153
Fig 4




**Warwick-James
UI-2-5154
Fig 1**



**Warwick-James
UI-2-5155
Fig 2**



**Warwick-James
UI-2-5156
Fig 3**



**London Hospital
UI-2-5157
Fig 1**




**London Hospital
UI-2-5158
Fig 2**



**London Hospital
UI-2-5159
Fig 3**



**Ott
UI-2-5160**



Pott



Pott
UI-2-5162
Fig 1



Pott
UI-2-5163
Fig 2



UI-2-5164
Fig 1



UI-2-5165
Fig 2



Mead
UI-2-5166
MD 1



Mead
UI-2-5167
MD-5 A




Mead
UI-2-5168
MD-6 A



Mead
UI-2-5169
MD-7 A



Mead
UI-2-5170
MD-8



Lecluse
UI-2-5171



**Chompret
UI-2-5172
Fig 1**



**Chompret
UI-2-5173
Fig 2**



**Chompret
UI-2-5174
Fig 3**



**Chompret
UI-2-5175
Fig 4**



**Apical
UI-2-5176**



UI-2-5177
Fig 1



UI-2-5178
Fig 2



UI-2-5179
Fig 3



Heidbrink
UI-2-5180
Fig 1



Heidbrink
UI-2-5181
Fig 2



Heidbrink
UI-2-5182
Fig 3



Winter
UI-2-5183
Fig 1L



Winter
UI-2-5184
Fig 1R



Winter
UI-2-5185
Fig 2L



Winter
UI-2-5186
Fig 2R



Winter



Winter
UI-2-5187
Fig 11 L



Winter
UI-2-5188
Fig 11 R



Winter
UI-2-5189
Fig 12 L



Winter
UI-2-5190
Fig 12 R



Winter
UI-2-5191
Fig 14 L



Winter
UI-2-5192
Fig 14 R



**Barry
UI-2-5193
Fig 1**



**Barry
UI-2-5194
Fig 2**



**Barry
UI-2-5195
Fig 3**



**Barry
UI-2-5196
Fig 4**



Barry



**Barry
UI-2-5197
Fig 11**



**Barry
UI-2-5198
Fig 12**



**Barry
UI-2-5199
Fig 17**



**Barry
UI-2-5200
Fig 18**



Cryer
UI-2-5201
Fig 26



Seldin
UI-2-5202
Fig 27



Cryer
UI-2-5203
Fig 28



Cryer
UI-2-5204
Fig 29



Heidbrink
UI-2-5205
Fig 22



Heidbrink
UI-2-5206
Fig 23



UI-2-5207

Root Elevators



Kopp
UI-2-5208
Fig 31



Kopp
UI-2-5209
Fig 32



Kopp
UI-2-5210
Fig 33



Kopp
UI-2-5211
Fig 34



Kopp
UI-2-5212
Fig 35



Kopp
UI-2-5213
Fig 44



UI-2-5214



Heidbrink
UI-2-5215
Fig 10



Winter
UI-2-5216
Fig 11



Winter
UI-2-5217
Fig 12



Barry
UI-2-5218
Fig 17



Barry
UI-2-5219
Fig 18



UI-2-5220
Fig 20



UI-2-5221

Root Elevators according to Dr. Sader,
with interchangeable working tips



UI-2-5222



UI-2-5223





UI-2-5224


Page 51-55

- Filling Instruments











UI-2-5268
Heidemann
Fig 1



UI-2-5269
Heidemann
Fig 2



UI-2-5270
Heidemann
Fig 3



UI-2-5271



UI-2-5272
Fig 146




UI-2-5274
Fig 148




UI-2-5273
Fig 147




UI-2-5275
Fig 149



UI-2-5276
Fig 158



UI-2-5277
Fig 159



UI-2-5278
Woodson
Fig 160



UI-2-5279
Gracey



UI-2-5280
Gracey
Fig 1



UI-2-5281
Gracey
Fig 2



UI-2-5282
Gracey
Fig 3



UI-2-5283
Gracey
Fig 4



UI-2-5284
Gracey
Fig 5



UI-2-5285
Gracey
Fig 6



UI-2-5286



UI-2-5287
Fig 18



UI-2-5288
Fig 19



UI-2-5289
Fig 20



UI-2-5290
Fig 22



UI-2-5291
Fig 23



UI-2-5292
Fig 25



UI-2-5293
Fig 26



UI-2-5294
Fig 27

Page 56–59

- Dental Probes



With Point Set In



UI-2-5295



UI-2-5296



UI-2-5297



UI-2-5298



UI-2-5299



UI-2-5300



UI-2-5301



UI-2-5302



UI-2-5303



UI-2-5304



UI-2-5305




UI-2-5306



UI-2-5307
Woodson

Polishing Instruments



Periodontal Pocket Probes



- New marking process.
- Laser system marking.
- No more paint peeling off.
- Permanent marking.
- Antiglare, low depth marking.
- Maximum tip resistance against breaking.



UI-2-5322
Fig CP-C




UI-2-5323
Fig CP-11



UI-2-5324
Fig CP-12



UI-2-5325
Fig CP-12S
O.M.S



Page 60-67

- Periodontia Instruments, Mallets



Files

Chisels



UI-3-5101
Sugarman
Fig. 1S/2S



UI-3-5102
Sugarman
Fig. 3S/4S



UI-3-5103
Schluger
Fig. 9/10



UI-3-5104
Buck
Fig. 11/12



UI-3-5105
Ochsenbein
Fig. 3




UI-3-5106
Ochsenbein
Fig. 4



UI-3-5107
Rhodes
Back - Action
Fig. 36/37

Files



UI-3-5108
Sugarman
Fig. 1S/2S

UI-3-5109
Sugarman
Fig. 3S/4S

UI-3-5110
Schluger
Fig. 9/10



UI-3-5111
Buck
Fig. 11/12

Chisels



UI-3-5112
Ochsenbein
Fig. 3



UI-3-5113
Ochsenbein
Fig. 4



UI-3-5114
Rhodes
Back - Action
Fig. 36/37

Gingivectomy knives



UI-3-5115
Goldman-Fox
Fig. 7



UI-3-5116
Goldman-Fox
Fig. 8



UI-3-5117
Goldman-Fox
Fig. 9



UI-3-5118
Goldman-Fox
Fig. 11



UI-3-5119
Buck
Fig. 3/4



UI-3-5120
Buck
Fig. 5/6

Curettes and Scalers



UI-3-5121
Curette
Fig. HU 1/2



UI-3-5122
Columbia
Fig. 2R/2L



UI-3-5123
Columbia
Fig. 4R/4L



UI-3-5124
Columbia
Fig. 13/14



UI-3-5125
Fig. 411/412



UI-3-5126
Fig. 413/414



UI-3-5127
Zerfing-chisel

Gingivectomy knives



UI-3-5128
Goldman - Fox
Fig. 7



UI-3-5129
Goldman - Fox
Fig. 8



UI-3-5130
Goldman - Fox
Fig. 9



UI-3-5131
Goldman - Fox
Fig. 11



UI-3-5132
Buck
Fig. 3/4



UI-3-5133
Buck
Fig. 5/6

Curettes and Scalers



UI-3-5134
Curette
Fig. HU 1/2



UI-3-5135
Columbia
Fig. 2R/2L



UI-3-5136
Columbia
Fig. 4R/4L



UI-3-5137
Columbia
Fig. 13/14



UI-3-5138
Fig. 413/412



UI-3-5139
Fig. 411/414



UI-3-5140
Zerfing-chisel

Bone Chisels and Gouges



UI-3-5141
Ochsenbein
Fig. 1

UI-3-5142
Ochsenbein
Fig. 2

UI-3-5143
Gardner
Fig. 1

UI-3-5144
Gardner
Fig. 2

UI-3-5145
Gardner
Fig. 5

UI-3-5146
Gardner
Fig. 5R


UI-3-5147
Gardner
Fig. 6

UI-3-5148
Wakefield
Fig. 1W


UI-3-5149
Wakefield
Fig. 2W

UI-3-5150
Buckley
Fig. 3

Bone Chisels and Gouges




Martelli chirurgici Mallet



Page 68-74

▪ Currettes Scalers





UI-3-5170
Gracey
Fig. 1/2



UI-3-5171
Gracey
Fig. 3/4



UI-3-5172
Gracey
Fig. 5/6



UI-3-5173
Gracey
Fig. 7/8



UI-3-5174
Gracey
Fig. 9/10



UI-3-5175
Gracey
Fig. 11/12



UI-3-5176
Gracey
Fig. 13/14




UI-3-5177
Gracey
Mesial Distal
Fig. 11/14



UI-3-5178
Gracey
Mesial Distal
Fig. 12/13

Gracey Curettes AF Type



UI-3-5179
Gracey
AF/type
Fig. 1/2



UI-3-5180
Gracey
AF/type
Fig. 3/4



UI-3-5181
Gracey
AF/type
Fig. 5/6



UI-3-5182
Gracey
AF/type
Fig. 7/8



UI-3-5183
Gracey
AF/type
Fig. 11/12



UI-3-5184
Gracey
AF/type
Fig. 13/14



UI-3-5185
Gracey
AF/type
Fig. 15/16



UI-3-5186
Gracey
AF/type



UI-3-5187
Gracey
AF/type
Fig. 1/2



UI-3-5188
Gracey
AF/type
Fig. 3/4



UI-3-5189
Gracey
AF/type
Fig. 5/6



UI-3-5190
Gracey
AF/type
Fig. 7/8



UI-3-5191
Gracey
AF/type
Fig. 11/12



UI-3-5192
Gracey
AF/type
Fig. 13/14



UI-3-5193
Gracey
AF/type
Fig. 15/16



UI-3-5194
Gracey
MF/type
Fig. 1/2



UI-3-5195
Gracey
MF/type
Fig. 3/4



UI-3-5196
Gracey
MF/type
Fig. 5/6



UI-3-5197
Gracey
MF/type
Fig. 7/8




UI-3-5198
Gracey
MF/type
Fig. 11/12



UI-3-5199
Gracey
MF/type
Fig. 13/14



UI-3-5200
Gracey
MF/type
Fig. 15/16



UI-3-5201
Gracey
MF/type
Fig. 1/2



UI-3-5202
Gracey
MF/type
Fig. 3/4



UI-3-5203
Gracey
MF/type
Fig. 5/6



UI-3-5204
Gracey
MF/type
Fig. 7/8



UI-3-5205
Gracey
MF/type
Fig. 11/12



UI-3-5206
Gracey
MF/type
Fig. 13/14



UI-3-5207
Gracey
MF/type
Fig. 15/16



UI-3-5208
Goldman - Fax
Fig. 1




UI-3-5209
Goldman - Fax
Fig. 2



UI-3-5212
Goldman - Fax
Fig. 21



UI-3-5210
Goldman - Fax
Fig. 3




UI-3-5213
Goldman - Fax
Fig. 5




UI-3-5211
Goldman - Fax
Fig. 4



UI-3-5214
Mc Call
Fig. 13S/14S



UI-3-5215
Younger-Good
Fig. 7/8





UI-3-5228
Hemingway
Fig. 1



UI-3-5229
Hemingway
Fig. 2



UI-3-5230
Hemingway
Fig. 3



UI-3-5231
Lucas
Fig. 85



UI-3-5232
Lucas
Fig. 86




UI-3-5233
Lucas
Fig. 87


Page
75-77

Modelling Instruments & Measuring Instruments





UI-4-5101 Le Cron UI-4-5102 Le Cron UI-4-5103 Le Cron UI-4-5104 MSY UI-4-5105 MSY UI-4-5106 Le Cron UI-4-5107 Hylin UI-4-5108 Hylin UI-4-5109 Hylin UI-4-5110 Zahle UI-4-5111 Zahle UI-4-5112 Zahle



UI-4-5113 Gritman UI-4-5114 Beale UI-4-5115 Roach UI-4-5116 Vehe UI-4-5117 P.K. Thomas Fig. 1 UI-4-5118 P.K. Thomas Fig. 2 UI-4-5119 P.K. Thomas Fig. 3 UI-4-5120 P.K. Thomas Fig. 4 UI-4-5121 P.K. Thomas Fig. 5



UI-4-5122
Boley




UI-4-5123
Zurich
130 mm



UI-4-5124
Vernier
12 cm



UI-4-5125
Vernier



UI-4-5126
Caroll



UI-4-5127
Castroviejo's
Scleral Marker



UI-4-5128
Iwanson



UI-4-5129
Albee



UI-4-5130
18 cm
7 1/8"



UI-4-5131
18 cm
7 1/8"



UI-4-5132
17 cm
6 3/4"
Gritman



UI-4-5133
17 cm
6 3/4"
Medicon



UI-4-5134
22 cm
8 3/4"



UI-4-5135
Antkoviak



UI-4-5136
Ermert



UI-4-5137
Ermert



Page 78–79

▪ Wax Knives/plaster



UI-4-5138
Fahnenstock
Fig. 1
12,5cm



UI-4-5139
Fahnenstock
Fig. 2
17cm



UI-4-5140
Lessmann
Fig. 1
12,5cm



UI-4-5141
Lessmann
Fig. 2
17cm



UI-4-5142
Con scodellino
Fig. 1
12,5cm



UI-4-5143
Con scodellino
Fig. 2
17cm



UI-4-5144
Fig, 1
17 cm



UI-4-5145
Gritman
Fig, 2
17 cm



UI-4-5146
Fig, 3
16 cm



UI-4-5147
Fig, 4
16 cm

A collection of dental instruments including pliers, a scaler, a probe, a mirror, and a pair of scissors, arranged vertically against a white background. A semi-transparent grey box is overlaid on the lower half of the image.

Page 80–85

Diagnostic, Explorer, Endodontic Instruments

Explorers



UI-4-5148
Fig 23



UI-4-5149
Fig 6



UI-4-5150
Fig 8



UI-4-5151
Fig 9



UI-4-5152
Fig 17




UI-4-5153
Fig 3
Rein



UI-4-5154
Fig 23/17



UI-4-5155
Fig 17



UI-4-5156
Fig.1
Spreader

UI-4-5157
Fig.2
Spreader

UI-4-5158
Fig.3
Spreader

UI-4-5159
Fig.4
Spreader


UI-4-5110
Fig.5
Spreader

UI-4-5111
Fig.6
Spreader

UI-4-5112
Fig.7
Spreader

Double-Ended Root Canal Explorers Plugger/Spreaders

Ø 0.4 Ø 0.4 Ø 0.6



Spreaders

Ø 0.6 Ø 0.5




Periodontal Pocket Probes

3,5 6 9 12

1 2 3 5 7 8 9 10 12 14



Root-canal pluggers




Endodontic set




UI-4-5132
Measuring Plate

Root-canal pluggers

Optimum Flexibility
Very Fine Points



Root-canal explorers



UI-4-5143
Fig.6/1

UI-4-5144
Fig.6/17

UI-4-5145
Fig.6/23

Root-canal excavators



UI-4-5146
Fig.31L

UI-4-5147
Fig.32L

UI-4-5148
Fig.33L

Page
86-87

Excavators

1 mm



UI-4-5149
Fig 5

1.2 mm



UI-4-5150
Fig 14

1 mm



UI-4-5151
Fig 17

1.2 mm



UI-4-5152
Fig 18

1.6 mm



UI-4-5153
Fig 19

2 mm



UI-4-5154
Fig 20

1 mm



UI-4-5155
ASH
Fig.243

1,4 mm



UI-4-5156
ASH
Fig.G2

0,8 mm



UI-4-5157
Fig.1

0,8 mm



UI-4-5158
Fig.2

1 mm



UI-4-5159
Fig.11

1 mm



UI-4-5160
Fig.12

1,2 mm



UI-4-5161
Fig.21

1,2 mm



UI-4-5162
Fig.22

1,6 mm



UI-4-5163
Fig.31

1,6 mm



UI-4-5164
Fig.32

2 mm




UI-4-5165
Fig.41

2 mm




UI-4-5166
Fig.42

A background image showing a dentist in a white coat and mask examining a patient's teeth. The patient is reclined, and the dentist is using a dental mirror. A semi-transparent grey box is overlaid on the left side of the image, containing the page number and title.

Page 88–89

Handles For Mouth Mirrors



Handles for Mouth Mirrors and Mouth Mirrors



Page 90–93

Spatulas





UI-5-5126




UI-5-5127



UI-5-5128



UI-5-5129



UI-5-5130




UI-5-5131





UI-5-5132



UI-5-5133







Page 94–100

- Needle Holders



UI-5-5150
Mayo-Hegar
16 cm



UI-5-5151
Mayo-Hegar TC
16 cm



UI-5-5152
Mayo-Hegar
18 cm



UI-5-5153
Mayo-Hegar TC
18 cm



UI-5-5154
MiUI-Ryder
18 cm



UI-5-5155
MiUI-Ryder
15 cm



UI-5-5156
Derf
12 cm



UI-5-5157
Hegar-Baumgarther
14 cm



UI-5-5158
Crile-Wood
15 cm



UI-5-5159
Crile-Wood TC
15 cm



UI-5-5160
Adson
17.5 cm



UI-5-5161
Olsen-Hegar
14 cm



UI-5-5162
Baumgarther
14.5 cm



UI-5-5163
Gillies
16 cm



UI-5-5164
Gillies
16 cm



UI-5-5165
Hegar-Baumgarther
14 cm



UI-5-5166
Toennis
18 cm



UI-5-5167
Toennis
18 cm



UI-5-5168
Lichtenberg
17 cm



UI-5-5169
Lichtenberg
17 cm



UI-5-5170
Eiselberg-Mathieu
19.5 cm



UI-5-5171
Boynton
12 cm



UI-5-5172
Arruga
16 cm



UI-5-5173
Boynton
12 cm



UI-5-5174
Arruga
16 cm



UI-5-5175
Castroviejo
14 cm



UI-5-5176
Castroviejo
14 cm



UI-5-5177
Castroviejo
14/60 cm



UI-5-5178
Castroviejo
18 cm



UI-5-5179
Castroviejo
18 cm



UI-5-5180
Castroviejo
14 cm



UI-5-5181
Castroviejo
14 cm



UI-5-5182
Castroviejo
14/60 cm



UI-5-5183
Castroviejo
18 cm



UI-5-5184
Castroviejo
18 cm



UI-5-5185
Mathieu
14 cm



UI-5-5186
Mathieu
14 cm



UI-5-5187
Mathieu
14 cm



UI-5-5188
Mathieu
14 cm



UI-5-5189
Mathieu
17 cm



UI-5-5190
Mathieu
17 cm



UI-5-5191
Mathieu
14 cm



UI-5-5192
Micro-Barraquer
18 cm



The tips have tungsten carbide inserts for better handling.



UI-5-5193
Mathieu
17 cm



UI-5-5194
Mathieu
17 cm



UI-5-5195
Micro-Barraquer
14 cm



UI-5-5196
Castroviejo
13 cm



UI-5-5197
Castroviejo-Barraquer
13 cm



UI-5-5198
Castroviejo
13 cm



UI-5-5199
Castroviejo
13 cm



UI-5-5200
18.03mm, cm



UI-5-5201
18.04mm, cm



UI-5-5202
18.05mm, cm



UI-5-5206
Fig. 6FR
Fraizer
15.02mm, cm













UI-5-5203
Fig. 12FR
Fraizer
15.04mm, cm



UI-5-5204
Fig. 10FR
Fraizer
15.03, 3mm, cm



UI-5-5205
Fig. 8FR
Fraizer
15.02, 7mm, cm



Page 101–109

- Syringes and Scissors



UI-5-5222

Water Syringe complete with
1 cannula (Lure-Lock)



UI-5-5223

Water Syringe complete with
1 cannula (Lure-Lock)



UI-5-5224

Metal water syringe, self-filling,
with 1 cannula (Lure-Lock)



UI-5-5225

Sparrel Barrel with
plunger suitable for
syringes



UI-5-5226

Cannula, only, with
Luer-Lock connector
suitable for syringes



UI-5-5227
Aspirating Syringe



UI-5-5228
Position Gauge



UI-5-5229
Breech Loading Syringe



UI-5-5230
Aspirating & Breech Syringe



UI-5-5231
Band Adaptor



UI-5-5232
Endo Finger Ruler
With Adjustable Slide



UI-5-5233
Aspirating Syringe




UI-5-5234



UI-5-5235



UI-5-5236



UI-5-5239
Tip Only, short



UI-5-5240
Tip Only, long



UI-5-5241



UI-5-5242

UI-5-5237 1,8 ml
UI-5-5238 2,2 ml
Cartridge Syringe complete
With 2 tips and without needles



UI-5-5243



UI-5-5244
Cork screw tip



UI-5-5245
Cork screw tip



UI-5-5246
Folding type, cork
screw tip



UI-5-5247
Folding type, cork
screw tip



UI-5-5248



UI-5-5249
Cork screw tip



UI-5-5250
Spare cap with
metric thread



UI-5-5251
With Whitworth
thread




UI-5-5252
Autoclavable Plastic Protection sleeve



UI-5-5253
Brass syringe with stainless steel moving parts




UI-5-5254
Brass syringe with stainless
steel moving parts



UI-5-5255
11.5 cm
UI-5-5256
13 cm
UI-5-5257
14 cm




Standard
Sharp / Sharp



UI-5-5261
11.5 cm
UI-5-5262
13 cm
UI-5-5263
14 cm




Standard
Blunt / Sharp



UI-5-5267
11.5 cm
UI-5-5268
13 cm
UI-5-5269
14 cm



Standard
Blunt / Blunt




UI-5-5273

UI-5-5274



Mayo
14.5 cm


Suture Scissors | Tijeras para Hilos | Fadenscheren



UI-5-5275



Metzenbaum-Delicate
14.5 cm




UI-5-5276 **UI-5-5277**


UI-5-5278



Metzenbaum
14.5 cm




UI-5-5279
Littauer
14 cm



Spencer
UI-5-5280 9 cm
UI-5-5281 13 cm




Bandage Scissors Lister
UI-5-5282 9 cm
UI-5-5283 11 cm
UI-5-5284 14 cm
UI-5-5285 18 cm



UI-5-5286 UI-5-5287 UI-5-5288 UI-5-5289




Iris
11.5 cm



UI-5-5290 UI-5-5291 UI-5-5292



Noyes / Castroviejo
12 cm




UI-5-52903
La Grange
11.5 cm



UI-5-5294
Universal
18 cm



UI-5-5295
Kelly
16 cm



UI-5-5296


UI-5-5297



Kelly
16 cm



UI-5-5298
Goldman-Fox
13 cm



UI-5-5299



UI-5-5300



Goldman-Fox
13 cm

Scissors, One Blade Serrated


- One blade is with fine serrations to avoid tissue or ligature slippage.
- Manufactured from high quality stainless steel to provide longer-lasting cutting edge.
- Single black ring for easy identification.




Super-Cut

- Cutting Blade horned to "razor" sharpness with opposing blade finished with fine serrations to avoid tissue slippage.
- Extremely smooth cutting with less trauma.
- Manufactured from high quality stainless steel to provide longer-lasting cutting edge.
- Single gold plated ring for easy identification.







Iris
11.5 cm




Metzenbaum
14.5 cm




Kelly
16 cm



UI-5-5322
La Grange
11.5 cm




Goldman-Fox
13 cm



Strabismus
11.5 cm


Micro Scissors



UI-5-5327
15 cm

UI-5-5328
15 cm

Micro Scissor
Stainless Steel



UI-5-5329
18 cm

UI-5-5330
18 cm

Micro Scissor
Stainless Steel




UI-5-5331
18 cm

Micro Scissor
Stainless Steel
for ligature, one blade serrated

Micro Scissors with Titanium Stainless Steel Tips




UI-5-5332
18 cm



UI-5-5333
18 cm

UI-5-5334
18 cm



UI-5-5335
15.5 cm

UI-5-5336
15.5 cm



Page 110-114

- Dental Tweezers





UI-5-5337
College
15 cm



smooth



UI-5-5338
College
15 cm



serrated



UI-5-5339
College
15 cm



smooth



UI-5-5340
College
15 cm



serrated



UI-5-5341
College
15 cm



smooth



UI-5-5342
College
15 cm



serrated



UI-5-5343
College
15 cm



smooth



UI-5-5344
College
15 cm



serrated

Dental Tweezers



UI-5-5345
Meriam
16 cm



smooth



UI-5-5346
Meriam
16 cm



serrated



UI-5-5347
Meriam
16 cm



smooth



UI-5-5348
Meriam
16 cm



serrated



UI-5-5349
Meriam
16 cm



smooth



UI-5-5350
Meriam
16 cm



serrated



UI-5-5351
Meriam
16 cm




smooth



UI-5-5352
Meriam
16 cm



serrated




UI-5-5353
College
14.5 cm, 5³/₄"

UI-5-5354
Flagg
16 cm, 6¹/₄"

UI-5-5355
Flagg
16 cm, 6¹/₄"

UI-5-5356
Landgraf
15 cm, 6"

UI-5-5357
15,5 cm, 6¹/₈"



UI-5-5358
15.5 cm, 6¹/₈"

UI-5-5359
London-College
15 cm, 6"

UI-5-5360
London-College
15 cm, 6"

UI-5-5361
London-College
15 cm, 6"

UI-5-5362
London-College
15 cm, 6"



UI-5-5363
CK1
15 cm, 6"

UI-5-5364
CK2
15 cm, 6"

UI-5-5365
GF2
16 cm, 6 1/4"

UI-5-5366
Gf1
16 cm, 6 1/4"

UI-5-5367
Kopp
15.5 cm, 6 1/8"



UI-5-5368
Kopp
15.5 cm, 6 1/8"

UI-5-5369
15 cm, 6"

UI-5-5370
16.5 cm, 6 1/2"

UI-5-5371
15.5 cm, 6 1/8"

UI-5-5372
15.5 cm, 6 1/8"



UI-5-5373
London-College
15 cm, 6"



UI-5-5374
London-College
15 cm, 6"



con cierre



UI-5-5375
Meriam
16 cm, 6 1/4"




serrated



smooth



UI-5-5376
Allen
18 cm, 7 1/8"



UI-5-5377
Albrecht
16 cm, 6 1/4"



UI-5-5378
Albrecht
16 cm, 6 1/4"



UI-5-5379
16 cm, 6 1/4"



UI-5-5380
16 cm, 6 1/4"




UI-5-5381
16 cm, 6 1/4"

Page 115-117

Rubber Dam Clamps



Rubber Dam Clamp Forceps



CLAMPS (PLIERS)





UI-5-5392
Winged
Labial Clamp



UI-5-5393
Winged
Edged Labial Clamp



UI-5-5394
Winged
Cervical Clamp



UI-5-5395
Winged
Labial Clamp

Premolars



UI-5-5396
Winged



UI-5-5397
Winged



UI-5-5398
Wingless

For Small Teeth, Elongated or With Narrow Collars



UI-5-5399
Winged



UI-5-5400
Winged



UI-5-5401
Wingless

For General Use With Adjustable Jaws According To Gingival Line Shapes



UI-5-5402
Winged



UI-5-5403
Edged, Winged



UI-5-5404
Wingless

For Premolars, With Elongated Collar



UI-5-5405
Winged



UI-5-5406
Edged, Winged



UI-5-5407
Wingless

For Premolars, With Elongated Collar



UI-5-5408
Wingless



UI-5-5409
Wingless



UI-5-5410
Winged

UI-5-5411
Wingless

For Small Lower Molars



UI-5-5412
Winged

UI-5-5413
Wingless

For Small Lower Molars With Round clamps



UI-5-5414
Winged

UI-5-5415
Wingless

For Square Lower Molars



UI-5-5416
Winged

UI-5-5417
Wingless

For Round Upper Molars



UI-5-5418
Winged

UI-5-5419
Wingless

For Molar Roots



UI-5-5420
Winged

UI-5-5421
Wingless

For Partially Broken Molars



UI-5-5422
Winged

UI-5-5423
Edged Winged

UI-5-5424
Wingless

For Round Molars, Partially Broken, Undersized or Inclined



UI-5-5425
Winged

UI-5-5426
Wingless



UI-5-5427
Wingless

UI-5-5428
Wingless

UI-5-5429
Winged

UI-5-5430
Winged

Page 118–122

- Matrix Retainer
- Matrix Bands





UI-6-5101
Ivory
63mm



UI-6-5102
Tofflemire-Universal



UI-6-5103
Tofflemire-Junior



UI-6-5104
Tofflemire-Senior



UI-6-5105
Ivory 8
5 mm



UI-6-5106
Ivory 8
6 mm



UI-6-5107
Ivory 8
7 mm



UI-6-5108
Ivory 8 N



UI-6-5109
Didner 8 P



UI-6-5111



UI-6-5112



UI-6-5113



UI-6-5114

1750-5 mm x 100 cm
1750-6 mm x 100 cm
1750-7 mm x 100 cm



UI-6-5110

For Molars and Premolars



UI-6-5115




UI-6-5116



UI-6-5117




UI-6-5118
Molar
Fig 7



UI-6-5119
Molar
Fig 8



UI-6-5120
Molar
Fig 9



UI-6-5121
Molar
Fig 10



UI-6-5122
Tofflemire



UI-6-5123
Tofflemire



UI-6-5124
Tofflemire



UI-6-5125
Bicuspid



UI-6-5126
Bicuspid



UI-6-5127
Bicuspid



UI-6-5128
Bicuspid



UI-6-5129
Bicuspid



UI-6-5130
Bicuspid



UI-6-5134



UI-6-5131
Molar



UI-6-5132
Molar





UI-6-5133
Molars



Page 120-122

Crown Removers
Scissors
Forceps
Amalgum Guns
Carries





Universal "HM"
12 cm, 4³/₄"



UI-6-5144



UI-6-5145



UI-6-5146
Beebee
10,5 cm, 4¹/₈"



UI-6-5147
Mershon



UI-6-5148
Morell



UI-6-5149
Morell



UI-6-5150
Morell



UI-6-5151
Morell



UI-6-5152
Treyman



UI-6-5153
Automatic Crown Remover
With 2 Points



UI-6-5154



UI-6-5155



UI-6-5156
Medium-Large
20 cm

UI-6-5157
Medium-Jumbo
20 cm

UI-6-5158
Large-Jumbo
20 cm

UI-6-5159
Jumbo-Jumbo
20 cm

UI-6-5160
Medium
15.5 cm

UI-6-5161
Large
15.5 cm

UI-6-5162
Jumbo
15.5 cm



Ø 2 mm
UI-6-5163
Straight
15 cm

Ø 2 mm
UI-6-5164
15 cm

UI-6-5165
15 cm
90°



UI-6-5166
Amalgam Well
Non Skid Base
Ø 3.5 mm

Page 123

- Copper Ring Removing Plierspaper Articulator Forceps

Copper Ring Removing Pliers, Paper Articulator Forceps

Copper ring removing pliers



UI-6-5167
Bohm
11,5 cm



UI-6-5168
Bohm
14,5 cm



UI-6-5169
Furrer
16 cm



UI-6-5170
Furrer
16 cm



UI-6-5171
Miller
15 cm



UI-6-5172
Miller
15 cm



UI-6-5173
Kelly
14 cm




UI-6-5174
Halstead-Mosquito
12 cm



UI-6-5175
Micro-Mosquito
10 cm


Page 124-127

▪ Rongeurs







UI-6-5176
Friedmann
14 cm, 5 1/2"



UI-6-5177
14 cm, 5 1/2"




UI-6-5178
Mead
14 cm, 5 1/2"



UI-6-5179
Mead
14 cm, 5 1/2"



UI-6-5180
Mead
17 cm, 6 3/4"



UI-6-5181
Schell
17 cm, 6 3/4"



UI-6-5182
Luer
15 cm



UI-6-5183
Friedman
14 cm



UI-6-5184
Mini Friedman
12 cm



UI-6-5185
Goldman-Fax
10 cm



UI-6-5186
Blumenthal
15 cm



UI-6-5187
Blumenthal 45
14,5 cm




UI-6-5188
Blumenthal 90
14 cm



UI-6-5189
Bohler
15 cm




UI-6-5190
Kerrison
15 cm




UI-6-5191
Blumenthal
16 cm, 6 1/4"




UI-6-5192
Blumenthal
16 cm, 6 1/4"



UI-6-5193
Blumenthal
16 cm, 6 1/4"



UI-6-5194
Kleinert-Kutz
14 cm, 5 1/2"




UI-6-5195
Ruskin-Mini
15 cm



UI-6-5196
Beyer
18 cm, 7 1/8"



UI-6-5197
Zaufal-Jansen
18 cm, 7 1/8"




UI-6-5198
Niro
No. 5 S, 14 cm, 5 1/2"



UI-6-5199
Cleveland
No. 5 S, 17 cm, 6 3/4"



UI-6-5200
No. 18, 15 cm, 6"



UI-6-5201
Cleveland
16 cm, 6 1/4"



UI-6-5202
Gardner
16 cm, 6 1/4"




UI-6-5203
Liston-mini
15 cm, 6"

Page 128–129


- Bone Elevators, Bone Files and Chisels





UI-6-5204
Williger
16 cm, 6 1/4"

UI-6-5205
Henahan
17 cm, 6 3/4"




UI-6-5206
Freer
19 cm, 7 1/2"

UI-6-5207
Bennet
17 cm, 6 3/4"

UI-6-5208
Seldin
19 cm, 7 1/2"

UI-6-5209
Seldin
19,5 cm, 7 3/4"

Bone Files and Chisels



Page 130–151

- Pliers For Orthodontics
- Mouth Gags
- Tongue Depressors
- Lip and Cheek Retractors
- Implantology Instrumetns



Pinze per formare loops
Loops forming pliers



UI-7-5101
Adams

Adams pliers for bending
wires up to 0,7mm. (.28").



UI-7-5102
Adams

Adams pliers for bending
wires up to 0,7mm. (.28").



UI-7-5103
Adams

Tweed loop forming pliers.
For forming loop in wires up
to 0,56 mm. (.22").



UI-7-5104
Adams

Aderer three jaws bending
pliers. For wires up
to 0,7 mm. (.28").



UI-7-5105
Adams

Aderer three jaws bending
pliers. For wires upto 0,7 mm.
(.28").

Ligature Cutter Pliers



UI-7-5106

Micro-miniature light wire cutter. It cuts soft ligature wires. For use in caso of limited access. Caution: maximum cutting capacity: soft ligature wire up to 0,30 mm. (.012").



UI-7-5107

Wire cutter. For ligature wires and pins. It is strong enough to bend archwire ends after distal cutting. Maximum cutting capacity up to 0,56 x 0,7 mm. (.022" x .028").



UI-7-5108

Wire cutter. For ligature wires and pins. It is strong enough to bend archwire ends after distal cutting. Maximum cutting capacity: up to 0,9 mm. (.036").



UI-7-5109

Wire cutter. For ligature wires and pins. It is strong enough to bend archwire ends after distal cutting. Maximum cutting capacity: up to 1,5 mm. (.044").

Pinze per formare laboratorio
Laboratory pliers



UI-7-5110

Flat nose pliers with smooth beaks for bending wires up to 0,7 mm. (.028").



UI-7-5111

Bending wires up to 0,7 mm. (.028").



UI-7-5112

Bending wires up to 0,7 mm. (.028").



UI-7-5113

Bending wires up to 0,7 mm. (.028").




UI-7-5114
Nance

Nance loop closing pliers, for bending of stop loops and special for Adams' clasps. Wires up to 0,7 mm. (.028").



UI-7-5115
Nance

Universal pliers for bending wires up to 0,9 mm. (.035") or cutting up to 0,7 mm. (.028").



UI-7-5116

Stainless steel pliers tray standard dimensions



UI-7-5117
Universal

Universal pin roughening, cutting and bending pliers. For wires up to 0,9 mm. (.035"). Will cut wires up to 0,7 mm. (.028").

Utility and specialty pliers



UI-7-5118
Reynolds

Reynolds contouring pliers.



UI-7-5119
Reynolds

Lingual arch pliers used to form double and triple back bends in 0,7 mm. (.028") or 0,9 mm. (.035") wire for insertion in lingual sheaths.



UI-7-5120
Weingart

Weingart pliers with serrated points.



UI-7-5121
De La Rosa

De La Rosa arch forming pliers. Designed to form and contour arch wires on round or rectangular wire. The beak has grooves of 0,4 mm (.016"), 0,46 mm (.018") and 0,56 mm (.022").



UI-7-5122
De La Rosa

Arch forming pliers with a non grooved contouring surface the forming of precise arches.

Loop formare pliers



UI-7-5123
Waldsachs

Waldsachs universal pliers for bending wires up to 1,0 mm. (.039") or cutting up to 0,7 mm. (.028").



UI-7-5124

Round and Concave pliers for loop forming and bending of thin round and square wires up to 0,5 mm. (.020").



UI-7-5125
Tweed

Tweed loop Forming Pliers. For round or rectangular wires up to 0,56 x 0,7 mm. (.022" x .028"). Replaceable cylindrical beak with three sizes: 1,15 mm. (.045), 1,50 mm. (.059"), 1,90 mm. (.075"). Concave beak with slight parallel serrations, to keep the wire in proper position.



UI-7-5126
Tweed

Tweed Loop Forming Pliers. For making perfect Omega Loops with round or rectangular wires up to 0,56 x 0,7 mm. (.022" x .028"). The concave beak is exactly the same as pliers 5809-1



UI-7-5127
Tweed

Single diameter loop forming pliers. Wire up to 0,56 x 0,7 mm. (.022" x 0.28"). Replaceable concave beak has only one diameter section of 1,50 mm. (.059").

Distal end cutters



UI-7-5128

Distal end cutter -safety hold- a precision cutter designed for distal arch cutting in the patient's mouth. It cuts hard wires end holds the loose distal end preventing any damaging of oral tissue. Maximum cutting capacity rectangular wires up to 0,56 x 0,7 mm. (.022" x .020").



UI-7-5129

Micro- distal end cutter - safety hold-. Similar to 3227



UI-7-5130

Direct bonding bracket remover (wide). The removal grip firmly for easy incisal-gingival (vertical) use or, if needed, for mesial-distal (horizontal) use with the same ease.




UI-7-5131

Direct bond attachment remover (angled). For removing all bonded attachments. May be used mesial-distal-horizontal) or incisal-gingival (vertical) in anterior areas. Angled perfectly for easy access to difficult posterior areas.



UI-7-5132

Wire bending pliers. For wires up to 0,9 mm. (.036").



UI-7-5133

Bird Beak Angle pliers for bending round, square or rectangular wire up to 0,7 mm. (.028").



UI-7-5134

Double rounded jaw wire bending pliers. For wires up to 0,7 mm. (.028").



UI-7-5135

Bird beak pliers for bending round, square or rectangular wire up to 0,7 mm. (.028"). With cutter.



UI-7-5136

Double rounded jaw wire bending pliers. For wires up to 0,7 mm. (.028").


Bird beack pliers



UI-7-5137
USA Model
125 mm




UI-7-5138
Waldron
125 mm



UI-7-5139
How
140 mm



UI-7-5140
How
145 mm




UI-7-5141
Adams
125 mm



UI-7-5142
Angle
125 mm



UI-7-5143
young
130 mm



UI-7-5144
Nance
130 mm




UI-7-5145
Weingart
140 mm




UI-7-5146
Johnson
130 mm




UI-7-5147
Mc Kellops
155 mm




UI-7-5148
Universal
150 mm



UI-7-5149
Waldsachs
155 mm



UI-7-5150
Aderer
125 mm



UI-7-5151
Nance
130 mm

Pliers for Othodontics and Prosthetics



UI-7-5152
Langenbeck
140 mm



UI-7-5153
Langenbeck
130 mm




UI-7-5154
Langenbeck
130 mm



UI-7-5155
Peeso
125 mm



UI-7-5156
Goslee
125 mm




UI-7-5157
Fischer
140 mm




UI-7-5158
Andresen
130 mm



UI-7-5159
Schwarz
130 mm



UI-7-5160
MSY
130 mm



UI-7-5161
Langenbeck
125 mm



UI-7-5162



UI-7-5163



UI-7-5164



UI-7-5165



UI-7-5166




UI-7-5167




UI-7-5168




UI-7-5169




UI-7-5170




UI-7-5171




UI-7-5172




UI-7-5173




UI-7-5174




UI-7-5175




UI-7-5176




UI-7-5177




UI-7-5178




UI-7-5179




UI-7-5180




UI-7-5181




UI-7-5182




UI-7-5183



UI-7-5184



UI-7-5185



UI-7-5186



Molt
UI-7-5187 10 cm
UI-7-5188 11 cm
UI-7-5189 14 cm



UI-7-5190 16 cm
UI-7-5191 19 cm
 Roser-König



UI-7-5192
 Denhardt
 13 cm



UI-7-5193
 Fergusson - Ackland
 14.5 cm



UI-7-5194
 Heister
 Maunder Wood
 9.5 cm




UI-7-5195
 Obwegeser 11 cm
 Spreading 47 mm



UI-7-5196
 Heister
 13 cm



UI-7-5197 Small
UI-7-5198 Medium
UI-7-5199 Large
 Cheek Retractor



UI-7-5200
Hartmann
15 cm



UI-7-5201
Doyen Tongue Spatula
13.5 cm




UI-7-5202
14.5 cm




UI-7-5203
18 cm




UI-7-5204
Bruenings
19 cm




UI-7-5205
Tobold
15 cm




UI-7-5206
Tobold
16 cm



UI-7-5207
14 cm



UI-7-5208
15 cm





UI-7-5209
Obwegeser Chin
Retractor
16.5 cm



Obwegeser 16 cm
UI-7-5210 8 mm
UI-7-5211 10 mm



UI-7-5212 43 x 11 mm
UI-7-5213 57 x 11 mm
UI-7-5214 83 x 16 mm

35 x 11 mm



UI-7-5215
Obwegeser
Spina-nasalis retractor
22 cm



72 x 11 mm

UI-7-5216
Obwegeser
Ramus retractor for
separation of ramus & tissue
22 cm



UI-7-5217
Simplex
10.5 cm



30x19 mm

33x25 mm

UI-7-5218
Kilner
15 cm



UI-7-5219
Minnesota
14 cm



UI-7-5220
Cawood-Minnesota
15 cm



UI-7-5221
Austin



UI-7-5222
Hilger
Fits any simple stem mouth mirror handle



UI-7-5223
Hilger
with handle



UI-7-5224
Max
20 cm

Finger protectors




UI-7-5225 Small
UI-7-5226 Medium
UI-7-5227 Large

Langenbeck



UI-7-5228
Oringer
11.5 cm



UI-7-5229
Obwegeser Cheek Retractor
12cm



UI-7-5230
Columbia
14.5 cm



Sternberg


UI-7-5231 12 cm
UI-7-5232 14 cm
UI-7-5233 15 cm



Vestibulum
UI-7-5234




UI-7-5235
Max




Implants Retractors

UI-7-5236
Max

Retractors




Senn-Miller
16 cm




Max
16 cm




Max
16 cm



Kirkland
Gum & Tissue
Retractors



UI-7-5249
Fig. 27 K



UI-7-5250
Fig. 28 K




UI-7-5251
Collin Tongue holding forcep
17 cm



Young Tongue holding forcep
UI-7-5252 15 cm
UI-7-5253 17 cm



UI-7-5255
Farabeuf 12 cm
Pair



UI-7-5256
Parker-Langenbeck 21 cm
Pair

Osteotomes



15.5 cm

# 1	Diam.
0	Ø 1.60
7.5	Ø 1.80
10	Ø 2.00
13.5	Ø 2.40
15.5	Ø 3.60

1
Ø 1.6 mm



UI-7-5257
15.5 cm

# 2	Diam.
0	Ø 1.90
7.5	Ø 2.20
10	Ø 2.35
13.5	Ø 3.30
15.5	Ø 3.40
23.5	Ø 4.50

2
Ø 1.9 mm



UI-7-5258
15.5 cm

# 3	Diam.
0	Ø 2.80
7.5	Ø 3.20
10.5	Ø 3.35
13.5	Ø 3.60
21	Ø 4.40

3
Ø 2.8 mm



UI-7-5259
15.5 cm

# 4	Diam.
0	Ø 3.10
7.5	Ø 3.70
10.5	Ø 3.90
15.5	Ø 4.30

4
Ø 3.10 mm



UI-7-5260
15.5 cm

# 5	Diam.
0	Ø 3.90
10	Ø 5.00
13	Ø 5.20
20.5	Ø 6.30

5
Ø 3.90 mm



UI-7-5261
15.5 cm

# 6	Diam.
0	Ø 5.50
10	Ø 6.20
13	Ø 6.30

6
Ø 5.50 mm



UI-7-5262
15.5 cm

# 1	Diam.
0	Ø 1.50
5	Ø 1.90
8	Ø 2.20
10	Ø 2.40
12	Ø 2.60
14	Ø 2.80

1
Ø 1.5 mm



UI-7-5263
15.5 cm

# 2	Diam.
0	Ø 2.10
5	Ø 2.60
8	Ø 2.80
10	Ø 3.00
12	Ø 3.20
14	Ø 3.40

2
Ø 2.1 mm



UI-7-5264
15.5 cm

# 3	Diam.
0	Ø 2.60
5	Ø 3.00
8	Ø 3.30
10	Ø 3.50
12	Ø 3.70
14	Ø 3.90

3
Ø 2.6 mm



UI-7-5265
15.5 cm

# 4	Diam.
0	Ø 3.00
5	Ø 3.50
8	Ø 3.80
10	Ø 4.00
12	Ø 4.20
14	Ø 4.40

4
Ø 3.0 mm



UI-7-5266
15.5 cm

# 5	Diam.
0	Ø 3.50
5	Ø 3.90
8	Ø 4.20
10	Ø 4.40
12	Ø 4.60
14	Ø 4.80

5
Ø 3.5 mm



UI-7-5267
15.5 cm

Ø 2.8 mm



UI-7-5268
15.5 cm

Spelling Straight
Bone Pusher

Ø 3.5 mm




UI-7-5269
15.5 cm

Ø 4.2 mm








UI-7-5270
15.5 cm

Osteotomes



Convex Tips





PDL Elevators



UI-7-5345 Distal, 2.5 mm Micro Edge



UI-7-5346 Mesial, 2.5 mm Micro Edge



UI-7-5347 Distal, 4 mm Micro Edge



UI-7-5348 Mesial 4 mm Micro Edge



UI-7-5349 2.5 mm Micro Edge, Straight



UI-7-5350 2.5 mm Micro Edge, Curved



UI-7-5351 3.4 mm Straight

UI-7-5352 Set Of 7 Pieces



Instruments Tray



UI-7-5353 Distal, 2.5 mm Micro Edge



UI-7-5354 Mesial, 2.5 mm Micro Edge



UI-7-5355 Distal, 4 mm Micro Edge



UI-7-5356 Mesial 4 mm Micro Edge



UI-7-5357 2.5 mm Micro Edge, Straight



UI-7-5358 2.5 mm Micro Edge, Curved



UI-7-5359 3.4 mm Straight

UI-7-5360 Set Of 7 Pieces

Trephine Drills



UI-7-5361
Interior 2 mm
Exterior 3 mm



UI-7-5362
Interior 3 mm
Exterior 4 mm



UI-7-5363
Interior 4 mm
Exterior 5 mm



UI-7-5364
Interior 5 mm
Exterior 6 mm



UI-7-5365
Interior 6 mm
Exterior 7 mm



UI-7-5366
Interior 7 mm
Exterior 8 mm



UI-7-5367
Interior 8 mm
Exterior 9 mm



UI-7-5368
Interior 9 mm
Exterior 10 mm



UI-7-5369
Trephine Drills Stand
Complete with drills



UI-7-5370
Trephine Drills Stand
Complete with drills

Tissue Punches



UI-7-5371
2 mm

UI-7-5372
3 mm

UI-7-5373
4 mm

UI-7-5374
5 mm

UI-7-5375
6 mm

UI-7-5376
7 mm

UI-7-5377
8 mm



UI-7-5378
3 mm

UI-7-5379
4 mm

UI-7-5380
5 mm

UI-7-5381
6 mm

Keyes

Keyes



UI-7-5382
Bone Collector Chisel Straight



UI-7-5383
Blade for Straight bone collector chisel



UI-7-5384
Bone Collector Chisel Curved



UI-7-5385
Blade for Curved bone collector chisel



UI-7-5386
Ratchet Screw driver handle
Hexagonal head



UI-7-5387
Bone Morselizer with Titanium Grinding



UI-7-5388
Bone Mixing Basin
Dia 35 mm



UI-7-5389
Manual Easy Bone Mill



UI-7-5390
Bone Craker Titanium



Page 152-155

- Impression Trays

All Impression Trays Upper & Lower Type
are available in S, M, L & XL Sizes



UI-7-5391
"HYDRO-LOCK"
stainless steel
water-cooled
impression trays
with retention rim.



UI-7-5392
"HYDRO-LOCK"
stainless steel
water-cooled impression
trays with retention rim.



UI-7-5393
"NEW SUPER"
stainless steel
impression trays
with retention rim.



UI-7-5394
"NEW SUPER" stainless
steel impression trays
with retention rim.



UI-7-5395
"ANATOMIC" solid
regular stainless
steel impression
trays without
retention rim.



UI-7-5396
"ANATOMIC" solid regular
stainless steel impression
trays without retention rim.



UI-7-5397
Perforated-total denture
stainless steel
impression trays with rim.



UI-7-5398
Perforated-total denture
stainless steel
impression trays with rim.

All Impression Trays Upper & Lower Type
are available in S, M, L & XL Sizes



UI-7-5399
"PERMA-LOCK" solid regular stainless
steel impression trays with retention rim.



UI-7-5400
"PERMA-LOCK" solid regular stainless
steel impression trays with retention rim.



UI-7-5401



UI-7-5402



UI-7-5403



UI-7-5404



UI-7-5405



UI-7-5406



UI-7-5407
Set of six

PACKINGS



UI-7-5408
Set of Eight



UI-7-5409
Solid stainless steel
anterior depressed
impression trays with
retention rim.



UI-7-5410



UI-7-5411



UI-7-5412



UI-7-5413



UI-7-5414



UI-7-5415



UI-7-5416
Solid stainless steel
anterior depressed
impression trays with
retention rim.

All Impression Trays Upper & Lower Type
are available in S, M, L & XL Sizes



UI-7-5417
Fig. 0 per bambini - for children



UI-7-5418
For narrow toothed upper jaws



UI-7-5419
Fig. 0 per bambini - for children



UI-7-5420
For partially toothed upper jaws



UI-7-5421
For partially toothed lower jaws



UI-7-5422
For edentulous upper jaws



UI-7-5423
For edentulous lower jaws

All Impression Trays Upper & Lower Type
are available in S, M, L & XL Sizes



UI-7-5424
Solid total denture stainless steel impression
tray with rim.



UI-7-5425
Solid total denture stainless steel impression
tray with rim.



UI-7-5426
"PERMA-LOCK" perforated regular stainless
steel impression trays with retention rim.



UI-7-5427
"PERMA-LOCK" perforated regular stainless
steel impression trays with retention rim.



UI-7-5428
Perforated stainless steel anterior depressed
impression trays with retention rim.



UI-7-5429



UI-7-5430



UI-7-5431



UI-7-5432



UI-7-5433



UI-7-5434



UI-7-5435
Perforated stainless steel anterior depressed
impression trays with retention rim.



UI-7-5436
Perforated stainless
steel impression trays
without retention rim
for complete denture.



UI-7-5437
Perforated stainless steel
impression
trays without retention
rim for complete denture.

Page 156–157

▪ Arriculators



UI-8-5101



UI-8-5102



UI-8-5103



UI-8-5104




UI-8-5105



UI-8-5106

Articulators



UI-8-5107



UI-8-5108



UI-8-5109




UI-8-5110

Page 158-160

- Dental Surgical Hollowares





UI-8-5111
ENDODONTIC BOX



UI-8-5112
STERILIZING BOX
Available in dif. sizes
17x17x3 cm
20x10x4 cm
25x12x6 cm
30x14x6 cm



UI-8-5113
ENDO TRAY BOX



UI-8-5114
COTTON DISPENSER



UI-8-5115
ALCOHOL LAMP

Sterilizing Drums



UI-8-5116



UI-8-5117



UI-8-5118

Sterilizing Boxes



UI-8-5119



UI-8-5120



UI-8-5121

Sterilizing Boxes



UI-8-5122



UI-8-5123



UI-8-5124



UI-8-5125



UI-8-5126



UI-8-5127

FORCEP AND THERMOMETER JARS

COTTON DISPENSER



UI-8-5128
IRRIGATOR(Douch can)



UI-8-5129
MEASURING JUG



UI-8-5130
KIDNEY TRAY



UI-8-5131
FORCEP JAR



UI-8-5132
GALLIPOT WITH & WITHOUT LID



UI-8-5133
SPITTEN MUG



UI-8-5134
SPIRIT LAMP



UI-8-5135
Perforated Tray



UI-8-5136
Perforated Base



UI-8-5137
Solid Base



UI-8-5138
Instruments Tray



UI-8-5139
Instruments Tray



UI-8-5140
**Instruments Tray
(Perforated)**



UI-8-5141
**Instruments Tray Shallow
(Mayo)**



UI-8-5142
Scalers Tray



UI-8-5143
Cathetar Tray



UI-101



UI-102



UI-103



UI-104



UI-105



UI-106



UI-107



UI-108



UI-109



UI-110



UI-111



UI-112



UI-113



UI-115



UI-114



UI-116




UI-117

Alphabetical Index

A			G		
Adams	125, 133	Cleveland	122	Goslee	134
Aderer	133	College	105, 107	Gouges	68
Adson	96	Columbia	64, 66, 74	Gracey Curettes	69
Ainsworth	110	Columbia	79	Gracey Mesial Distal	69
Albee	77	Cone-Socket	90	Gracey	56-69-72
Albrecht	109	Cork Screw Tip	104	Gritman	76-78-80
Alcohol Lamp	142	Cotton Dispenser	142, 143		
Alginate	93-94	Coupland	37	H	
Allen	109	Crane-kaplan	74	Halstead-Mosquito	118
Amalgam Carriers	117	Crile-Wood	96	Hegar-Baumgartner	95-96
Amalgam Guns	117	Crown Removers	115-116	Heidbrink	45-48-50
Amalgam Well	117	Cryer	35, 37, 48	Heidemann	55
Andresen	134	Curette	64, 66	Hemingway	75
Angle	133			Henahan	123
Antkoviak	78	D		How	133
Apical	36, 44	De La Rosa	129	Hylin	40, 76
Arruga	97	Dental Probes	57		
Articulators	140	Dental Syringes	103-104	I	
Ash	53, 88	Derf	95	Impression Trays	136-139
Aspirating Syringes	102	Didner	113	Intra-Ligamental Syringes	104
				Irrigator	143
B		E		Ivory	110, 113
Band Adaptor	102	Eliselberg-Mathieu	97	Iwanson	77
Barry	47, 50	Eliott	53		
Baumgarthner	96	Endo Finger Ruler	102	J	
Beale	76	Endo Tray Box	142	Jaquette	74
Beebee	116	Endodontic Box	142	Johnson	133
Bein	33, 40	Ermert	78		
Bennet	53, 123	Excavators	87-88	K	
Beyer	121	Explorers	81	Kelly	118
Bicuspids	114	Extracting Forceps	01-31	Kerrison	120
Birdbeak Pliers	132			Kidney Tray	143
Blumenthal	120, 121	F		Kilian	124
Bohler	120	Fahnenstock	79	Kilian-Claus	124
Bohm	118	Filling Instruments	52-56	Klein	07
Boley	77	Fischer	134	Kleinert-Kutz	121
Bone Chisels	68	Flagg	107	Kopp	39, 49, 108
Boynton	97	Flohr	38		
Breech Loading Syringe	102	Forcepps Jar	143, 144	L	
Brewer	110	Fox	58	Landgraf	107
Buck	61-63, 65	Fraizer	100	Langenbeck	134
Buckley	67, 124	Freer	123	Le Cron	76
Burnishers	58	Friedman	119-120	Lecluse	43
		Furrer	119	Lessmann	79
C				Lichtenberg	97
Caroll	77	G		Listion-mini	122
Castroviejo	77, 98, 100	Gallipot	144	London Hospital	41
Cathetar Tray	144	Gardner	67, 122	London-College	107, 109
Cement Spatulas	91, 92	GF1	108	Lucas	75
Chompret	44	GF2	108	Luer	120
Ck1	108	Gillies	96	Luks	84
CK2	108	Goldman Fox	60, 63, 65, 73, 120		

M			W		
Machtou	83	Rhodium	89	Wakefield	
Mallets	68	Roach	76	Waldron	67
Mathieu	99, 100	Root Elevators	32--51	Waldsachs	133
Matrix Bands	113, 114	Root Splinter	31	Warwick-James	130, 133
Matrix Retainers	113	Root-Canal Excavators	86	Water Syringe	41
Mayo-Hegar	95	Root-Canal Explorers	86	Weingart	101
Mc Call	73	Root-Canal Plungers	84-85	White	129, 133
Mc Kellops	133	Rubber Dam Clamps	111-112	Williams	52, 60, 83
Mead	43, 58, 119	Rubber Dam Frame	110	Williger	123
Measuring Jug	143	Ruskin-Mini	121	Winter	46, 50
Measuring Plate	84			Woodson	54, 55
Medicon	78	S			
Meriam	106, 109	Scalers	58	Y	
Merrit A	60	Schell	119	Young	133
Mershon	110, 116	Schilder	85	Younger-Good	73
Micro-Barraquer	99-100	Schluger	61, 62	Z	
Micro-Mosquito	118	Schwarz	134	Zahle	76
Miller	118, 124	Seldin	34, 48, 123	Zaufal-Jansen	121
MiUI-Ryder	95	Spirit Lamp	144	Zerfing-Chisel	64, 66
Molar	114	Spitten Mug	144	Zuricher	77
Morell	116	Spreader	82-83		
Mouth Mirror Handles	89, 90	Sterilizing Boxes	142, 143		
Mouth Mirrors	89, 90	Sterilizing Drums	142		
MSY	76, 134	Sugarman	61, 162		
N			T		
Nabers	83	Taylor	74		
Nance	128, 133	Toennis	97		
Niro	122	Tofflemire	114		
O			Tofflemire-Junior	113	
Obwegeser	124		Tofflemire-Senior	113	
Ochsenbein	61, 62, 67		Tofflemire-Universal	113	
Olsen-Heger	96		Towner Jaquette	74	
Oms	83		Towner	74	
Ott	41		Trays	144	
P			Treymann	116	
P.K. Thomas	76		Tweed	130	
P.K. Thomas	68	U			
Partsch	134	Univarsal "HM"	116		
Pocket Probes	58-60	Univarsal	128, 133		
Plaster Knives	80	University of Michigan	60		
Plaster Spatulas	93-94	USA Model	133		
Position Gauge	102				
Pott	42	V			
R			Vehe	76	
Rein	81		Vernier	77	
Relax	05		Vertical Condenser	83	
Reynolds	129				
Rhodes	61, 62				



UNZILA INDUSTRY

ALL TYPES OF DENTAL INSTRUMENTS

contact us

Mobile No : 0092 3086403561

Web Site : www.unzilaindustry.com

E-mail : unzilaindustry@gmail.com

Address : Ugoki road near sobedar street Sialkot

